

GETINGE SMART DISTRIBUTION TROLLEYS
SMART MOVES

GETINGE

GETINGE GROUP



Always with you

DISTRIBUTION MADE SMARTER FOR DIFFERENT NEEDS AND DIFFERENT USERS

Getinge SMART Distribution Trolleys are designed for the effective, ergonomic transportation of baskets and/or containers. They are used in, to and from packing areas in sterile zones, and for transportation to soiled zones throughout the hospital. hospital - simply for the easy transfer of goods. The trolleys' development has been underpinned by ergonomical thinking, best usability and quality design.

DISTRIBUTION TROLLEY HARD COVER (DTHC)

In the DTHC instrument containers and baskets are grouped together with utensils, linen, disposables and anesthetic material to be transported to an internal OR-department, or external hospital/facility storage. The hard cover protects the sterile goods from dust and damage during transportation.

- Simple and safe transportation
- Ergonomic
- Can be automatic washed and disinfected

Capacity

8 pcs SPRI basket (585x395x195)

15 pcs DIN (600x300x155)

DISTRIBUTION TROLLEY OPEN (DTO)

The DTO is a mobile storage unit for modular baskets and containers. The open design makes the viewing and picking of individual packs very easy. With shelves and baskets easily secured for safe transportation.

- Easy to load and unload
- Saves storage space
- Silent and easy maneuverable

Capacity

10 pcs SPRI basket Low

5 pcs SPRI basket High

10 Guides applicable for shelves



DISTRIBUTION TROLLEY HARD COVER



The doors can be opened in 260° and remain in open position, which makes loading and unloading very convenient. Door handles are prepared for safety lock.



As an option, the trolley can be equipped with additional shelves for extra storage space. Each shelf is delivered with 2 dividers which are enabling the goods to remain in position during transport.



To facilitate the manoeuvring of the DTHC it is equipped with 4 wheels, of which 2 are lockable and 2 has directional lock. Optional pin and tow bar is available for transport after truck.



DISTRIBUTION TROLLEY OPEN



The flexible placement of support guides and shelves maximizes the utilization of the trolleys.



The loading is available from both sides for easier loading and unloading. You can also pull the shelf/basket out half way for easier access.



The high quality disinfectable wheels are even when fully loaded easy to move and maneuver, and for safety two of the wheels are lockable.



COMPLETE SOLUTIONS FOR INFECTION CONTROL

Getinge is the world's leading provider of solutions for effective cleaning, disinfection and sterilization in the healthcare and life science sectors. We are dedicated to helping our customers provide better care at a lower cost. We do this by offering well thought through and customized solutions. This means that we are with our customers all the way from architectural planning and education to traceability and support – with complete solutions, long-term commitment and global presence. Getinge – Always with you.

Getinge Infection Control AB

P O Box 69
SE-310 44 Getinge
SWEDEN
Phone: +46 10 335 00 00
info@getinge.com

GETINGE

GETINGE GROUP

www.getinge.com

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. **ArjoHuntleigh** focuses on patient mobility and wound management solutions. **GETINGE** provides solutions for infection control within healthcare and contamination prevention within life sciences. **MAQUET** specializes in solutions, therapies and products for surgical interventions and intensive care.