



Sterility Assurance Products  
World Class and Globally Trusted



## Sterilization Control Tubes

### Sterilization Control Tubes

#### Type 1 Black spot

Ref: 7301

Fluid Sterilization - 121°C/15 mins - Class 4

#### Type 2 Yellow spot

Ref: 7302

Steam Sterilization\* - 134°C/3 mins - Class 4

#### Type 3 Green spot

Ref: 7303

Dry Heat Sterilization - 160°C/60 mins - Class 6

#### Type 4 Blue spot

Ref: 7304

Dry Heat Sterilization - 180°C/12 mins - Class 6

#### Type 5 White spot

Ref: 7305

Dry Heat Sterilization - 160°C/120 mins - 180°C/35 mins - Class 6


Conform to ISO 11140-1<sup>2</sup>

Effective and immediate visual check on sterilization conditions


- Integrated 2 criteria control of temperature and time
- Ideal for fluid and dry heat sterilizers

Select the Control Tube to match your sterilization parameters (time/temperature) and distribute them throughout the load prior to sterilization.


Time to go Green (mins). Type 1 Black Spot Tubes




Time to go Green (mins). Type 2 Yellow Spot Tubes




Time to go Green (Mins). Type 3 Green Spot Tubes



Time to go Green (Mins). Type 4 Blue Spot Tubes



Time to go Green (Mins). Type 5 White Spot Tubes



### Colour Tubes Scale



Unused



Unsafe



Effective Cycle

Sterilization

fluid, dry heat and steam sterilization